

# SMF CoffeePRO

SYSTEM PART:	#80-6100CP
REPLACEMENT CARTRIDGE:	#108-010CP
RATED CAPACITY (GALLONS):	CHLORINE *: 50,000 CHLORAMINE: 15,000 SCALE INHIBITION: 15,000
FLOW RATE:	5 GPM
MIN./MAX. PRESSURE:	40– 125 PSI
MIN./MAX. TEMP.:	40°– 100° F
INLET/OUTLET SIZE:	1/2" BARB STANDARD 17 MODULAR FITTINGS AVAILABLE
SHIPPING WEIGHT:	17 LBS.
DIMENSIONS (W/D/H):	14/9/ 19"

## COFFEE & ESPRESSO APPLICATIONS

FOR: SCALE INHIBITION, CHLORINE,  
CHLORAMINE, SEDIMENT, THMS,  
ORGANICS, TANNINS, BAD TASTE AND  
ODOR REDUCTION

VOLUME: MEDIUM



### SPECIFIED TO:

- INHIBIT SCALE BUILD-UP
- REDUCE PITTING OF HEATING ELEMENTS IN BREWING EQUIPMENT
- IMPROVE TASTE OF BREWED BEVERAGES
- PROVIDE SAFE, CLEAN INGREDIENT WATER
- REDUCE CHLORINE, CHLORAMINE, BAD TASTES AND ODORS, TOTAL ORGANIC COMPOUNDS, TANNINS AND TRIHALOMETHANES

### FILTER BENEFITS:

- ONE FILTER CAN DO IT ALL: DIRT REDUCTION, CHEMICAL REDUCTION, CHLORAMINES REDUCTION AND SCALE INHIBITION! NO PRE-FILTER REQUIRED!
- THE FIRST JUMBO SANITARY FILTER ON THE MARKET!
- NYLON CONSTRUCTION TO WITHSTAND HIGH PRESSURE!
- FILTER PERFORMANCE GUARANTEE!
- 40 TIMES MORE DIRT HOLDING CAPACITY THAN CONVENTIONAL CARBON BLOCK FILTERS AND OTHER POLYPROPYLENE SEDIMENT FILTERS!
- 30 TIMES THE CHLORAMINE REDUCTION OF A CONVENTIONAL CARBON BLOCK!
- LOW PRESSURE DROP!
- MODULAR FILTER DESIGN FOR UP TO 50% LESS INSTALLATION TIME!
- QUICKLY ADD FITTINGS FOR EXTRA FILTER CAPACITY WITH NO TOOLS REQUIRED!

### TECHNOLOGY:

- ULTRA HIGH POROSITY CARBON WITH SILIPHOS

### APPLICATIONS:

- COFFEE BREWERS
- ESPRESSO & CAPPUCCINO MACHINES

\*Independently tested by an E.P.A. approved laboratory to ANSI/NSF Standard 42 specifications for the reduction of Chlorine, Chloramine, Taste and Odor.



Base system has been tested and certified by NSF International against ANSI/NSF Std. 42 Particulate Class I, 0.5-1 µm. Please visit [www.nsf.org](http://www.nsf.org) for specifications.

Do not use with water that is microbiologically unsafe or of unknown quality.

**selecto**

3980 LAKEFIELD COURT, SUWANEE, GA 30024 USA  
800-635-4017 • 678-475-0799 • [WWW.SELECTOINC.COM](http://WWW.SELECTOINC.COM)

© 2008 Selecto Scientific, Inc. • All Rights Reserved.

805-009 Rev. 03-27-08