

SMF Espresso 7500

SYSTEM PART:	#82-1600
REPLACEMENT CARTRIDGE:	#108-600
RATED CAPACITY (GALLONS) @ 5 GPM:	CHLORINE *: 10,000 CHLORAMINE *: 3,000 SCALE REDUCTION (GALLONS) : 7,500
FLOW RATE:	1 GPM
MIN./MAX. PRESSURE:	40– 125 PSI
MIN./MAX. TEMP.:	40°– 100° F
INLET/OUTLET SIZE:	1/2" BARB STANDARD 17 MODULAR FITTINGS AVAILABLE
SHIPPING WEIGHT:	19 LBS.
DIMENSIONS (W/D/H):	14/9/ 19"



COFFEE, TEA & ESPRESSO APPLICATIONS

FOR: SCALE REDUCTION, CHLORINE,
CHLORAMINE, SEDIMENT, THMS,
ORGANICS, TANNINS, BAD TASTE AND
ODOR REDUCTION

VOLUME: LOW-MEDIUM

SPECIFIED TO:

- INHIBIT SCALE BUILD-UP
- REDUCE PITTING OF HEATING ELEMENTS IN COFFEE/TEA BREWERS
- IMPROVE TASTE OF BREWED BEVERAGES
- PROVIDE SAFE, CLEAN INGREDIENT WATER
- REDUCE CHLORINE, CHLORAMINE, BAD TASTES AND ODORS, TOTAL ORGANIC COMPOUNDS, TANNINS AND TRIHALOMETHANES

FILTER BENEFITS:

- ONE FILTER CAN DO IT ALL: DIRT REDUCTION, CHEMICAL REDUCTION, CHLORAMINES REDUCTION AND SCALE REDUCTION! NO PRE-FILTER REQUIRED!
- THE FIRST JUMBO SANITARY FILTER ON THE MARKET!
- NYLON CONSTRUCTION TO WITHSTAND HIGH PRESSURE!
- FILTER PERFORMANCE GUARANTEE!
- 40 TIMES MORE DIRT HOLDING CAPACITY THAN CONVENTIONAL CARBON BLOCK FILTERS AND OTHER POLYPROPYLENE SEDIMENT FILTERS!
- 30 TIMES THE CHLORAMINE REDUCTION OF A CONVENTIONAL CARBON BLOCK!
- LOW PRESSURE DROP!
- MODULAR FILTER DESIGN FOR UP TO 50% LESS INSTALLATION TIME!
- QUICKLY ADD FITTINGS FOR EXTRA FILTER CAPACITY WITH NO TOOLS REQUIRED!

TECHNOLOGY:

- CERAMIC MATRIX WITH ULTRA HIGH POROSITY CARBON
- NANOCRYSTAL

APPLICATIONS:

- COFFEE BREWERS
- TEA BREWERS
- ESPRESSO & CAPPUCCINO MACHINES

*Independently tested by an E.P.A. approved laboratory to ANSI/NSF Standard 42 specifications for the reduction of Particulate Class II, Chlorine, Chloramine, Taste and Odor.



Base system has been tested and certified by NSF International against ANSI/NSF Std. 42. Please visit www.nsf.org for specifications.

Do not use with water that is microbiologically unsafe or of unknown quality.

selecto

3980 LAKEFIELD COURT, SUWANEE, GA 30024 USA
800-635-4017 • 678-475-0799 • WWW.SELECTOINC.COM

© 2007 Selecto Scientific, Inc. • All Rights Reserved.

805-010 Rev. 10-18-07